

DATA SHEET WORKFORCE - First Aid & Bleed Control Trauma Kit

SKU Number SS-6010125060 Date of Publication 09-01-2025

EAN13 Barcode 5060660858194 Commodity Code 30065000.00

Product Image







WORKFORCE - First Aid & Bleed Control Trauma Kit (complies with BS 8599-1:2019)

When severe haemorrhage control is needed, having a kit equipped with the right medical products is essential. Simple, user-friendly Bleed Control Kits are critical to the effectiveness of first responders.

Product Description

The **STEIN WORKFORCE** is a hybrid kit for First Aid and Bleed Control Trauma Kit. Combining the elements of a FRONTLINE Plus Kit with that of a Small Workplace First Aid Kit complying with BS 8599-1:2019; this kit is packed with high-quality medical supplies necessary to manage life-threatening haemorrhages, lacerations, burns, and minor injuries.

This kit allows trained First Aiders, Advanced Professionals, and First Responders to provide lifesaving care in critical situations whilst waiting for Emergency Services to arrive.

Contained in a durable weather resistant bag made from advanced materials, this kit is designed for quick access in emergencies. The elements of the First Aid and Bleed Control are separated into two separate internal bags distinguished by colour.



	Orange Bag:
	1 x Tactical Multi Bandage 1 x SAM SPLINT 9" Flat 1 x High Performance Emergency Thermal Blanket 3 x Security Seals
	Green Bag (Fist Aid):
	Small Workplace First Aid Kit conforming to BS 8599-1:2019
Contents	Red Bag (Bleed Control):
	1 x SAM XT Tourniquet 1 x SWAT-T Rescue Orange Tourniquet 1 x 6" FCP Military Emergency Bandage 1 x ChitoSAM Haemostatic Gauze 4' Z-Fold 1 x Sterile Saline Spray (100ml) 1 x High Performance Medical Shears 1 x Medical Marker Pen 2 x Disposable Medical Gloves 2 x De-Contamination Wipes 30cm x 20cm
Safety Standards	All medical products used in the STEIN Kits are CE Marked and conform with the relevant European Union harmonisation legislation.
Product Weight	• 2300g
Size	• 320cm x 270cm x 120cm (approx.)
Storage	 Keep stored and transported at ambient temperature (0°C to + 30°C) Do not store in direct/high heat or sunlight as this may affect the internal components of the kit
Lifetime	 Not all products have an expiry date, however most sterile items do. If a non-sterile product has no expiry date and still sealed in its original packaging, it's a question of judgement as to whether it is still suitable for use. Items that have exceeded their expiry date should be disposed and replaced as they may not be effective and can start to deteriorate, making them unreliable.
Cleaning	Clean the storage bag using warm soapy water and soft cloth
Disposal	Recycling is preferred where possible. Look out for the recycling symbol for material categories



Replacement Components	 Following the use of any product the kit should be replenished with items used so they are fully functional for future incidents Refer to the STEIN Website for replacement items www.steinworldwide.com
Additional Information	 In extreme cases of haemorrhage further products may be required over and above the items in these kits. This kit does not replace the mandatory Workplace First Aid Kits complying with BS 8599-1:2019
Training	 It is not feasible to address every potential scenario involving the use of this product. Users are strongly advised to obtain professional training from a fully qualified and competent instructor before use. Training Providers can be recommended upon request
Trademarks	STEIN® is a registered trademark